


# Skin training suture for small animal



With the various pre-set wounds, this pad offers the optimal training basis for beginners and advanced beginners.

**Product codes:**

Reference: VTQ-SUT-PAD-SA

**Product description:**

This three-layered suture pad is particularly soft and flexible, perfectly mimicking the characteristics of delicate dog and cat skin.

It is suitable for practicing incisions, sutures and surgical knot techniques under realistic conditions.

The visual and textural demarcation between skin and subcutaneous tissue additionally provides opportunity for intracutaneous suture practice.

Our suture pads are reusable, tear-resistant and retain their shape and texture even after repeated use.

Our practical pad holder (S) can be used to mount suture pads in a vertical or a horizontal position, as required.

**Product features:**

Dimensions: 15 x 16 cm

### Available models

