

# DECRO Rat nude jacket

protective Rat jacket (w/o emitter)



## Product codes:

Reference: JKT-RAT-NU

## Product description:

DECRO jackets are available in a nude version, featuring the same design as the instrumented jackets but without embedded sensors.

Nude jackets are notably suited :

- to protect catheters or head-fixed implants
- protect wounds from the animal,
- or to embed other type of emitter (activity trackers, etc) .

They can also be used for ECG monitoring by integrating DECRO ECG electrodes. In this configuration, they function as part of the DECRO solution, ensuring seamless compatibility with DECRO.

Size to be confirmed when placing the order

**Product features:**

**Available models**

