

# Cordless trimmer for rodent BIO-FIDELA



Robust cordless trimmer, offering excellent handling for rodent

**Product codes:**

Reference: BIO-FIDELA

**Product description:**

Cordless trimmer, equipped with a quiet and powerful motor, allowing very fine cuts on a wide variety of animals.

Super light and quiet, only 141g. Ideal for handling rats and mice.

Timing charge : 150 minutes

Delivered with its surgical cutting head **reference BIO-TCFIDELA**

**Product features:**

Width: 30 mm

Weight: 141 gr

Charging time: charge : 150 min

Cutting length: 0.4 mm  
Battery: 120 min

#### Available models

