

Autoclip Wound Closing System 9mm- Full Kit



Maximum holding force with minimal trauma

Product codes:

Reference: ACS-KIT

Product description:

The Autoclip® Applier works like **a staple gun to rapidly apply Autoclip closures**. The applier will hold 20 Autoclips, which are held firm and ready for use by rustless beryllium-copper springs. Autoclips provide maximum holding force with minimal trauma. There is **no need for assistance to operate the applier or remover**.

- When the wound is healed, clips are easily removed with the Autoclip Remover.
- Stainless steel Autoclips are non-pyrogenic and supplied non-sterile, but may be autoclaved.
- Ideal for use with laboratory mice and rats.

Product features:

Available models