

Mini Incubator



Compact and economically, yet offering a wide range of functionalities

Product codes:

Reference: I5110

Product description:

Compact and economically priced, offers features not typically found in a basic incubator.

Product features:

Weight: 8.3 kg/18.4 lb

Power supply: Controls : Analog

Temperature range: Temperature range

Dimensions - Inner: 9 x 7.9 x 7.9 in (23 x 20 x 20 cm)

Dimensions - Outer: 11.2 x 11 x 13.2 in (28.5 x 28 x 33.5 cm)

Volume: 9.2 liters/0.375 cu. ft. 6

Specification: Electrical : 120V~ or 230V~, 50/60 Hz

Available models