

## Mini LabRoller Rotator

Economical solution for tube mixing



### Product codes:

Reference: H5500

### Product description:

Simple, yet versatile, the Mini LabRoller Rotator is designed for mixing samples in both a horizontal and vertical plane. The unit accepts a broad assortment of tube sizes as well as micro plates.

Safe for cold room or incubator use, this mixer is ideal for a wide range of applications in areas such as biochemistry, molecular biology, clinical testing and histochemistry.

### Product features:

Dimensions: (WxDxH) : 8.4 x 4 x 5 in 21.3 x 10.2 x 12.6 cm

Working temperature: 4° to 40°C

Speed: 115V/230V : 24 / 20rpm

Specification: Motion/Speed : Rotating, rocking, tumbling

### Available models