


Intrathecal/Intercranial Catheter for Rats

Intrathecal Catheter



Product codes:

Reference: R-ITC

Product description:

Rat intrathecal catheter offers a means for the **administration of microliter amounts of drug solution to the lumbar** (or higher) region of the rat spinal cord.

This catheter system is ideal for use in **antinociceptive studies** involving rats.

Product features:

Material: MicroRenathane[®] Tubing

Other features: for needle type N-27/27-PP

Length: total : 140mm / tip : 90mm

Diameter: O.D. 0.25mm

Available models