


## Double Lumen Catheters



Get two access points with a single vessel puncture or cutdown ! Extra small models available

### Product codes:

Reference: CR-DL

### Product description:

#### 3 different models are available:

- **DL-AS-040:** This double lumen catheter has an outer diameter of only 0.040" and consists of an inner SUBL-160 (PTFE) catheter and an outer MRE040 (polyurethane) catheter. The standard length of this double lumen catheter is 6" with both lumens opening at the proximal tip. The distal ends are equipped with luer hubs. Supplied non-sterile.
- **CR-DL28:** Double lumen intracranial catheter system contains an 11.7cm 28g (0.025" OD distal portion/0.007" ID each lumen) polyurethane catheter, a bonded connection assembly flared to accept a 27ga needle, a 0.003" diameter Teflon®-coated stainless steel stylet to facilitate placement. Individually packaged. Supplied non-sterile.
- **CR-DL32:** Double lumen intracranial catheter system contains an 11.7cm 32ga (0.018" OD distal portion/0.005" ID each lumen) polyurethane catheter, a bonded connection assembly flared to accept a 32ga needle, a 0.003" diameter Teflon®-coated stainless steel stylet to facilitate placement. Individually packaged. Supplied non-sterile.

**Product features:****Available models**

DL Catheter Models: 11.7 cm 28g (0.025" O.D. distal portion/0.007" I.D. each lumen) (CR-DL28),  
11.7 cm 32g (0.018" O.D. distal portion /0.005" I.D. each lumen) (CR-DL32)

