


# Universal Restrainer



Experiments Made Easier - different models for rats and mice, with smoked or clear canvas

**Product codes:**

Reference: 500M-C

**Product description:**

The Universal Restrainer was designed to allow the handler to perform all sampling and injection procedures in a single immobilization.

- the slit in the abdominal ventral area of the rat allows : access to the bladder, gastrointestinal catheters, fistulas, intraperitoneal administration, cardiac punctures
- the dorsal slit allows : the insertion of catheters, subcutaneous injections, EKG and EEG
- the slit at the tail allows : for sampling and injection, heart puncture

Made of Stainless steel and plastic

**Product features:**

Weight: mice 20g to 60g

Weight: rat up to 350g

**Available models**

Color: Clear (500M-C), Smoked (500M-S), Clear (700R-C), Smoked (700R-S)

Animal: Mice 20g to 60g (500M-C), Mice 20g to 60g (500M-S), Rats (less than 350g) (700R-C),  
Rats (less than 350g) (700R-S)

