


Restrained Tail Vein Infusion Model



An easier access for both timed injections and repeated bolus administration

Product codes:

Reference: MTV-

Product description:**Product features:****Available models**

Size: 10cm of 28ga (0.36x0.18mm) connected to 17.8 cm of (0.84x0.36mm) with stylette and 23ga thin wall needle (MTV-1), 10cm of 28ga (0.36x0.18mm) connected to 50.8cm of (0.84x0.36mm) and 23ga thin wall needle, no stylette (MTV-2), 10cm of 28ga (0.36x0.18mm) connected to 50.8cm of (0.84x0.36mm) with stylette and 23ga thin wall needle (MTV-2S), 10cm of 28ga (0.36x0.18mm) connected to 80cm of (0.84x0.36mm) and 23ga thin wall needle, no stylette (MTV-3)