

BioLITE Intubation Illumination System



Provides better anatomical visualization and combines several tools into one reducing the intubation time making the process fast and easy

Product codes:

Reference: BIO-KIT

Product description:

If you thought intubating small rodents was too difficult, the BioLITE Intubation Illumination System is for you. This system provides **better anatomical visualization** and combines several tools into one reducing the intubation time making the intubation process **fast and easy**.

The BioLITE System consists of a family of products: **Fiber Optic Illuminator (FO-Light)**, **Fiber Optic Stylets (FOS)**, **Endotracheal Tubes (ETT)**, and **Lung Inflation Bulbs (LIB)** packaged in a variety of kits to meet the individual needs of a variety of research interests.

Product features:

Available models

BioLITE Intubation System: Mice 10+ g (includes 1 FOS-400, 1 FO-Light, 1 LIB-03 and 1 ETT-20-05) (BIO-MI-KIT), Small Rats 100-200 g (includes 1 FOS-600, 1 FO-Light, 1 LIB-03 and 1

ETT-16-05) (BIO-SR-KIT), Large Rats 200+ g (includes 1 FOS-800, 1 FO-Light, 1 LIB-03 and 1 ETT-16-05) (BIO-LR-KIT), Stylets Kit (includes 1 FOS-400, 1 FOS-600, 1 FOS-800, 1 FO-Light) (FOS-KIT), Complete kit (includes 1 FOS-Kit, 1 LIB-03, 1 ETT-16-05 and 1 ETT-20-05) (BIO-CO-KIT)

