

## Elisabethan collars



Designed to prevent animals from grooming and to protect surgical or treatment sites

**Product codes:**

Reference: EC-MR

**Product description:**

All collars are manufactured from a low-density polyethylene, for smaller animals a lightweight material is used.

The collars are closed with Velcro, which makes it possible to adjust them for each animal.

The neck openings are bound by a padded cotton jersey for comfort.

**Product features:****Available models**

Animal: Mice 70mm (EC-M1), Rats up to 450g - 117mm (EC-R1), Rats 450g+ - 130mm (EC-R2)