


# Polyurethane Rat Vascular Catheters



New Polyurethane (Micro-Renathane<sup>®</sup> and RenaPulse<sup>®</sup>) Rat Vascular Catheters

**Product codes:**

Reference: R-CAC-

## Product description:

- Catheters consist of a continuous length of MRE037 or RPT037 with RenaSil, suture anchors.
- Catheters are accessed with Cat # N-23 catheter needle and sealed with Cat #23-PP catheter plug.
- No adhesives; suture anchors grip the catheter tubing without constricting the catheter lumen and are movable.
- Catheter tip rounding available for MRE037 catheter fabrications.

## Product features:

## Available models

Composition: MRE-037 (R-CAC-M37), RPT-037 (R-CAC-R37), MRE-037 - Tip rounding (R-CAC-M37-R)