


Tethers



Stainless steel spring-coil rodent tethers

Product codes:

Reference: RT

Product description:

Tethers Constructed from lightweight stainless steel spring-coil, rodent tethers come complete with a crimped end for direct attachment to the rodent swivel (removing the need for potentially troublesome connection hardware). Jacket anchors comes ready sewn with silver woven fabric for added comfort and resistance to opportunistic infection. Rodent tethers are designed to protect the exteriorized catheter and can be used with rodent jackets.

Product features:**Available models**

Size: 12" (RT-12), 18" (RT-18)