

Small Rodents Ear Tag



Small Animal Ear Tag, especially suited for Research Lab Animals, like mice, rats and other rodents

Product codes:

Reference: EP-1005-1

Product description:

Extensive and proven identification technique, ear tags have the advantage of their **simplicity**, in their laying, but also in monitoring the identified animals.

Product features:

Material: Monel (Nickel-Copper Alloy)

Stamping: 3 digits consecutive stamped number will fit on one side - 6 digits is available but will extend around other side

Available models

Ear tags: Engraved from 1 to 1,000 (EP-1005-1), Engraved - Numbering to be defined (EP-CUSTOM), Custom numbering & lettering sequence (EP-CUSTOM-L), Custom Lasar printing – number sequence (EP-CUSTOM-LR), blank (EP-CUSTOM-B)