


# Surgery and MicroTattooing Platform



General micro-utility platform ideal for surgery, micro manipulations, ...

**Product codes:**

Reference: NTS01

**Product description:**

MicroTattooing Platform with Single Flexible LED Light can be used for surgery, micro manipulations, general micro-utility platform (i.e. organ inspections, dissections, etc.).

- Flex arm, non-obstructive rimless magnifying glass (main lens: 2.3 X; inset lens: 4.0 X)
- Non-skid feet for added stability.
- Universal input power supply (100-240 VAC, 50-60 Hz).
- Surgery Platform incorporates one flexible LED light which provides direct illumination, no shadows and no heat.
- Includes 1.5mL cup holders.
- Size: 6"x9" (15.24cm X 22.86Ccm)

**Product features:****Available models**

