

## GyroMini Nutating Mixer



An ideal solution for mixing samples in blood and centrifuge tubes in a three dimensional action

**Product codes:**

Reference: S0500

**Product description:**

The GyroMini is a smaller version of the popular GyroTwister. Its small size makes it ideal for mixing samples in blood and centrifuge tubes.

**Product features:**

Weight: 2.3 lb / 1.1 kg

Dimensions: (W x D x H) : 9.75 x 7.5 x 7 in 24.8 x 19.1 x 17.8 cm

Temperature range: Ambient : + 4° to 65°C

Speed: Range : 24rpm

Specification: Maximum load : 1.75 lb / 0.8 kg

**Available models**