


Pie Cage for X-ray irradiation



An easy to use and to clean solution for temporary transportation and storage

Product codes:

Reference: PC-

Product description:

Circular cage holds **up to 11 rodents**.

An individual chamber is available to insert the rodents in the cage. All pie-shape chambers are ventilated.

Cages are available for mouse, small and big rats.

Ref

Specie

Diameter (cm)	Height (cm)	Indiv. Chamber (cm)	MPC-1	Mice	21	5,9	5 (width) x 9 (length)
			MPC-2	Mice	idem	idem	with filters
			RPC-1	Small Rat	41	6,75	10 (width) x 18 (length)
			RPC-2	Large Rat	46	8	11,5 (width) x 20 (length)

Product features:

Available models

Application: MPC-1 for mice (MPC-1), MPC-2 for mice, with filters (MPC-2), RPC-1 for small rats (RPC-1), RPC-2 for large rats (RPC-2)

