

# Rectal thermometer for rodents with footswitch



BIOSEB R&D team Innovation: finally a thermometer that leaves your hands free when taking temperature on rats and mice! An ideal tool to deal with large cohorts of animals (rats and mice) in a semi-automatic way. The perfect instrument for pharmacology, physiology monitoring, fever evaluation, inflammation monitoring, and much more!

**Product codes:**

Reference: BIO-THERMO-CIS

**Product description:**

THERMO-CIS solution was developed with a view of convenience to allow the experimenter to deal with large cohorts of animals by sending the data in an Excel spreadsheet of your choice and keep your **hands free** when handling animals.

Their **intuitive** operation, allow the thermometer and dedicated software to be the **ideal solution for temperature measurement routine of a large number of animals.**

**Product features:**

Power supply: 220V-5V 200mA provided (12 volt battery rescue)

Units: °C, °F

Resolution: Grade II thermocouple +/- 0,5°C (certificate on request)

Other features: Backlight : YES

Measuring range: -20 - + 200°C

Connection: USB 2,3 cable up to 3 meters

Specification: PC : Laptop, tablet computer with 1 USB port

## Available models

