

# Starter kit RapID Tags

The **RapID Tag Starter Kit** contains:

- 200 RapID tags (10 trays - 5 colors)
- 40 RapID practice tags (1 tray)
- 1 RapID applicator
- 1 RapID removal tool

## Product codes:

Reference: BIO-RAPIDSTART



## Product description:

The **RapID Tag System** has been designed to provide **safe, humane, automated tracking and data collection** for small laboratory animals.

RapID Tag's unique design encapsulates a **hardened metal barcode in a tough, chew-resistant, medical-grade, advanced polymer body**. The virtually indestructible, scratch-resistant barcode contains an exclusive combination of alpha-numeric characters providing the lab technician with **millions of unique ID numbers**.

Multi-colored RapID Tags also provide **ultra-fast (secondary) visual identification** without disturbing caged animals. RapID Tags can be automatched to computer-generated, colored cage-card labels using an exclusive website application.

RapID Lab Product Guide Video from RapID Lab on Vimeo.

## Product features:

## Available models

