

GentleSharp System

A more humane blood sampling system



Product codes:

Reference: GS-101-2

Product description:

GentleSharp is designed for venipuncture and arterial puncture of superficial blood vessels, with collection of blood via a capillary tube, dried blood spot card, or other disposable.

Handpiece allows for attachment of all **Luer Slip and Luer Lock needles** commonly used in blood sampling.

- Easy clean exterior and ergonomic handgrip ensure ease of use.
- RoHS compliant.
- Powered via electrical outlet (240V AC).
- Manufactured in ISO 13485 certified (cGMP compliant facility)

Product features:

Available models

Cable: 61cm (2") (GS-101-2), 122 cm (4") (GS-101-4)

