

# Rodent Sling and Restraint Cover



The best choice to securely restrain rodents during many procedures

**Product codes:**

Reference: RSL

**Product description:**

Slings can be used to securely restrain rodents during many procedures. Choice of restraint in many sizes and range of adjustment make rodent sling a valuable asset to any laboratory setting.

Other small animal species may also be restrained using same frame with suits and slings available for mice and guinea pigs.

**Product features:****Available models**

Description: Rodent Sling - frame and cover (RSL-FS1), Rodent Sling Cover Only (RSL-FC01), Rat Sling Suit Small, 150-200g (RSL-FC1), Rat Sling Suit Medium, 250-350g (RSL-FC2), Rat Sling Suit Large, 350-500g (RSL-FC3), Rat Sling Suit X-Large, 500g+ (RSL-FC4), Mouse Sling Suit (RSL-C00), Mouse Sling Suit Large, 30g+ (RSL-C03)