

Compact stereotaxic frame for mice



Small dimensions - wide working area

Product codes:

Reference: BLAB-ST-8507

Product description:

This stereotaxic instrument is designed to use adapter block to replace the traditional U-frame assembly, and it holds the mouse or neonatal rat firmly and accurately. Ear bars are made of polyoxymethylene resin and can be fixed using 4 different points with 1mm resolution for scales. 2 reversible mouse ear bars, with atraumatic tips, are sold with this product: 18° and 45°.

This compact instrument provides easy access to the animal and offers an open working area for extra surgical instruments.

Product features:

Dimensions: 255x255 mm

Available models

Manipulators: Single (BLAB-ST-8507), Dual (BLAB-ST-68510), Single (BLAB-ST-68513), Dual (BLAB-ST-68516)

Frame type : Manual (BLAB-ST-8507), Manual (BLAB-ST-68510), Digital (BLAB-ST-68513),
Digital (BLAB-ST-68516)

