

# Deltaphase Board



This compact operating board takes advantage of the Deltaphase principle to maintain animal body temperature during surgical or experimental procedures.

**Product codes:**

Reference: 39OP

**Product description:**

This compact operating board takes advantage of the Deltaphase principle to maintain animal body temperature during surgical or experimental procedures.

**An stainless steel operating surface lies in contact with an activated Deltaphase® isothermal pad.**

The underside of the pad is insulated to provide a longer isothermal period. Molded wells are provided for sponges, instruments and anesthetic apparatus. The board is functional, durable, and easy to clean. Its low cost and safety make it appropriate for use in student laboratories.

The unit is supplied with two Deltaphase® pads.

**Product features:**

## Available models

