

## Rat restrainer for BP-2000



to be used on BP-2000 unit only

### Product codes:

Reference: BP-RH

### Product description:

Restrainer for rat for BP-2000

Size-adjustable in length, for an ideal handling of the rat

Feature a magnet-baseplate, for an accurate fixation on the heating platform

Black color and holes at the front side keep the animal quiet for the duration of the blood pressure test.

Easy to clean and urine-resistant

### Product features:

Material: anodized metal

Colors: black

**Available models**

Animal Size: Rat as large as 1 kg (BP-RH6), Rat 700-850 grams (BP-RH5), Rat 550-700 grams (BP-RH4), Rat 400-550 grams (BP-RH3), Rat under 250 grams (BP-RH1), Rat 250-400 grams (BP-RH2)

