

## Endcaps for BP-2000



Bag of 10 endcaps

**Product codes:**

Reference: BP-CAP

### **Product description:**

One endcap is used at each end of the cuff to secure the cuff elastic.

The endcap size depends on the diameter of the cuff : it is recommended to check the diameter of the metal cuff (1/4", 3/8", 5/8") before to proceed with an order.

One bag contains 10 endcaps

### **Product features:**

### **Available models**

Size: 5/8" for rat (BP-CAP-R), 3/8" for mice (BP-CAP-M-38), 1/4" for mice (BP-CAP-M-14)