

# Tissue forceps



Standard Pattern Tissue Forceps. Secure ribbed grip. Available for standard (standard SS) or micro-tissues (premium grade stainless steel and narrow fine micro tips).



## Product codes:

Reference: FCT

## Product description:

## Product features:

## Available models

Tissue forceps: FC1-10 Standard pattern, Serrated Handle, 1x2 Teeth, 5" (FC110), FC1-09 Premium, serrated handle, 1x2 dent, 4.5" (FC109), FC1-11 Standard, Serrated Handle, 1x2 Teeth, 4.5" (FC111), FC14-5 Micro-Tissue forceps, Delicate Tissue, Curved Tip, 4.25" (FC145), FC14-6 Micro-tissue forceps, Delicate Tissue, Narrow Fine Tip, 4.25" (FC-146), FC14-7 Micro-tissue forceps, Delicate Tissue, Super Fine Tip, 4.25" (FC-147)