

## Metzenbaum scissors



Offering a unique combination of narrow blunted points along with thin delicate shafts for delicate dissection and minimal damage to surrounding tissues. Available in Standard or SuperCut version, curved or straight.

The SuperCut models are manufactured with the highest grade surgical stainless steel. Tempered and hardened for a longer sharper life. The SuperCut blade has a specially designed razor sharp edge.

### Product codes:

Reference: SCM

### Product description:

### Product features:

### Available models

Metzenbaum Scissors: SC-M-126 Curved, Premium Grade, 7" (SC-M-126), SC-M-125 Curved, Premium Grade, 5" (SC-M-125), SC-M-182 Curved, Standard Pattern, 7" (SC-M-182), SC-M-180 Curved, Standard Pattern, 5.5" (SC-M-180), SCT-M-525 Supercut, Curved, 4.5", delicate (SCT-M-525), SCT-M-523 SuperCut, Straight, 4.5", delicate (SCT-M-523)