


Deltaphase Isothermal Pads



Safe, reusable heating pads - requiring no electricity neither cable. Maintain animals or cultures at 37° for hours: a very useful tool for maintaining animals' temperature (rats or mice) during surgical procedures or to speed up recovery processes. Ideal for NMR (Naked Mole Rats)!

Sold by unit

Product codes:

Reference: DPIP

Product description:

Based on fundamental thermodynamics, completely safe even with flammable gases, uses no electricity, no wires, no control box, and they are reusable.

The Deltaphase Isothermal Pad is a new approach to body temperature stabilization. It acts as a source of heat, but cannot overheat. Its constant temperature is based on a thermodynamic principle and needs no electricity or wires, no thermostat or feedback controller. It cannot generate electrical signals or cause interference. It cannot trigger an explosion, but it can maintain a small animal at near normal and constant temperature for several hours.

Product features:

Available models

Size: 6" x 4" (DP-SASB), 6" x 7,5" (DPIP)