

Disposable DETRUNCATION restrainers



The economical restrainer for rodents

Product codes:

Reference: DC-00

Product description:

Make injections and decapitation quicker and easier with DecapiCones. These tapered plastic film tubes provide quick and easy restraint of rats, mice, and other small animals. I.P. injections can be made directly through the film! DecapiCones restrain post-decapitation kicking and prevent personal contact with feces or urine.

Product features:**Available models**

Description: Mice (Includes: 2 Dispensers, 100 bags each Total 200 cones) (MDC-200), Rats (Includes: 4 Dispensers, 50 bags each Total 200 cones) (DC-200), For larger rats (500-800g) (Includes: 4 Dispensers, 30 bags, Total 120 cones) (DCL-120)