

# HL-1 infrared Heat Lamp



a sturdily IR lamp for heating subjects after surgery procedures

**Product codes:**

Reference: HL-1

**Product description:**

Infrared heat lamp with an **adjustable swing arm providing for easy adjustment of radiant heat to specimen.**

The lamp is supplied with :

- a 250W(infrared) heat lamp bulb
- a built in reflector
- A chrome bulb guard is also provided and easily snaps onto the lamp housing after the bulb is installed.
- Arm Length: 2 spring loaded arms, each 45,7cm lengths

**HL-1** : Standard lamp is provided with a **Table Clamp**: Standard "U" screw clamp expands to fit bench up to 5cm deep. .

**HL-1/BASE** : Infrared heat lamp with **tabletop base**.

**Product features:**

**Available models**

Application: Table Clamp "U" screw (HL-1-LAMP), with tabletop base (HL-1-BASE)

