

## Red Induction box

Reassuring color for rodents



**Product codes:**

Reference: AB-RED

### Product description:

Red Gas Anesthetizing Box **with perforated floor** for use with rodents.

- The 0,6cm thick plastic anesthetizing box allows for constant observation while anesthetizing - Help protect against accidental overdosing
- **2 sizes** : For use with rats and mice
- Easy to move machine tooled sliding top
- **Removable interior, perforated floor** permits the rodent to be elevated off the floor allowing for a cleaner environment
- Made of durable plastic, this box should last for years
- Not for use with liquid organic solvents

### Product features:

### Available models

AB-: S (25.4 cm L x 10.16 cm W x 10.16 cm H) (AB-RED1), L (25.4 cm L x 20.3 cm W x 20.3 cm H) (AB-RED2)

