

# Flexible Disposable Feeding Needle with Bite protector



box of 100 sterile, flexible PTFE Disposable Feeding Needle, with bite protector

**Product codes:**

Reference: BG-99

**Product description:**

Plastic (PTFE) animal feeding needles (AFN)

with **protective Stainless steel hub**, avoid needle damage due to biting during use

allow user to precisely dispense liquid to animal while **minimizing the chances of esophageal piercing**.

sterile, non-toxic, and non-pyrogenic, available in a variety of gauges and lengths.

**Product features:**

Units: Each box has 100 sterile needles

**Available models**

Size: 17G 3.3" 2.4mm (BG-9926), 18G 3" 2mm (BG-9924), 18G 2" 2mm (BG-9923), 18G 1,5" 2mm (BG-9928), 15G 3.9" 2.8mm (BG-9927), 15G 3.1" 2.8mm (BG-9925)

