

Gas filter canister



Pack of 6 canisters for the adsorption of residual anesthetic gases

Product codes:

Reference: R510-31-6

Product description:

These gas filter canister are made of activated carbon, a material with a porous structure, which captures waste anesthetic gases to reject a "clean" air in the room. It connects directly to your induction cage or anesthetic masks via a bellows. Its lifespan is controlled by weighing. Once the defined weight is reached, the filter is considered saturated and needs to be changed. These canisters are sold in pack of 6.

Product features:

Units: Pack of 6

Norms compliance: CE certification

Available models

Size: Small (net weight: 250g) x 6 (R510-31-S6), Large (net weight: 800g) x 6 (R510-31-6)