

# Economy Flat Bottomed Restrainers



Semicylindrical plastic restrainers available in different sizes

**Product codes:**

Reference: FB-

**Product description:**

Each unit has an open slot on the bottom (for intra-peritoneal injections), a hole on the top (for the neck), and a slotted tail piece (for intravenous injections). Front slots provide ample ventilation

**Product features:****Available models**

Animal: Mice 15 to 30g (1.5" D x 4") (FB-S), Rats 30 to 125g (2" D x 5") (FB-M), Rats 125 to 250g (2.5" D. x 6") (FB-M/L), Rats 250 to 500g (3.375" D. x 8.5") (FB-L)