

Anesthesia cone mask

Allows access to eyes, ears and top of head

**Product codes:**

Reference: 68683

Product description:

Those cone masks are made of transparent acrylic. The diaphragm of these masks allows adjustment around the muzzle of animals and ensures a fair fit. They also have an inlet tubing inside that allow the animal anesthesia exhaust gas to be absorbed by gas filter canister. This eliminates potential anesthetic gas leaks in the room.

Product features:

Material: Clear acrylic

Available models

Version: Cone mask with tubing for neonatal mouse (less than 15g) (68680), Cone mask with tubing for mouse or neonatal rat (15-40g) (68681), Cone mask with tubing for rat (200-350g) (68682), Cone mask for neonatal mouse (less than 15g) (68683), Cone mask for mouse or neonatal rat (15-40g) (68684), Cone mask for rat (200-350g) (68685)