

Active anesthetic masks for stereotaxic frames

Fit most stereotaxic frame



Product codes:

Reference: BLAB-ST-8665

Product description:

Those anesthetic masks are specifically designed to fit most stereotaxic frame. They are locked onto the adaptor with the cone opening toward the open "U". The animal's teeth is placed over the incisor bar.

This kind of mask is called "active" ; it can thus be used with our gas evacuation apparatus. It has separate (non-collinear) inlet and outlet on both sides.

Product features:

Available models

Animal: Mice or neonatal rats (less than 30g) (BLAB-ST-8663), Mice or neonatal rats (30-70g) (BLAB-ST-8665), Rats (less than 300g) (68666), Rats (more than 300g) (68667)