

# USB Rechargeable trimmer for animals BIO-LA9600VET



Cordless, handy, light, powerful trimmer with 2 batteries - Recommended for big animals (cat, dog...)

## Product codes:

Reference: BIO-LA9600VET

## Product description:

Quiet, heavy duty, motor-driven rechargeable trimmer equipped with super fine tooth blades for close trimming.

2 lithium batteries : operate for 4 hours (2x2h)

2 blades set stainless steel included :

- ref TC9600 33 teeth
- ref TC9600P 5 teeth

Charging cable with USB plug

Cutting lenght adjustement with 5 positions 1 -3 mm by 0,5mm

4 attachment combs : 3 mm - 6 mm - 9 mm - 12 mm

Charging light green or red for battery indication

To make trimming very small animals and narrow areas easier, we offer an optional precision cutting head part number **TC9600P**.

Special application : for rodent, we recommend the trimmer BIO-1584

#### **Product features:**

Width: 38 mm

Weight: 200gr

Cutting length: 0.25mm without comb

#### **Available models**