

Precision cutting head for trimmer LA9600VET



Allows precise and narrow-cut , compatible with trimmer LA9600VET

Product codes:

Reference: BIO-TC9600P

Product description:

Precision cutting head for trimmer LA9600VET

Allows very precise shearing of your rodents,

- for the ears,
- tighter fit
- ideal for surgical procedures

Stainless steel

Cutting blade 5 teeth

Special application : for rodent we recommend the trimmer BIO-1584

Product features:

Cutting length: 1 mm
Trim head width: 8 mm

Available models

