

## Tailveiner for rat



Easy to use rat restrainer specifically designed for tail vein injections and samplings

**Product codes:**

Reference: RTV

**Product description:**

The Tailveiner is the first easy-to-use rat restrainer specifically designed for tail vein injections and samplings

The base is fixed by 3 suction cups to the bench, which allows a good stability.

- RTV-170 : up to 120g ID : 4,45cm
- RTV-180 : up to 250g ID : 6,98cm
- RTV-190 : up to 500g ID : 8,25cm

Comes with a stopper to stop the animal moving.

Made of plastic, not autoclavable, wash only soap and warm water

The red tube is transparent and available with all items.

**Product features:**

Weight: up to 120g

Weight: up to 500g

Weight: up to 350g

### Available models

Color: Red (RTV-R-170), Red (RTV-R-180), Red (RTV-R-190), Clear (RTV-170), Clear (RTV-180), Clear (RTV-190)

Rat Size Tailveiner : 4,45 cm up to 120g (RTV-R-170), 6,98 cm up to 350g (RTV-R-180), 8,25 cm up to 500g (RTV-R-190), 4,45 cm up to 120g (RTV-170), 6,98 cm up to 350g (RTV-180), 8,25 cm up to 500g (RTV-190)

