

Anesthesia Adaptor

for a stereotaxic frame

**Product codes:**

Reference: 68078

Product description:

This adaptor holds stably the rodent head by mask and incisor bar. The connection between the head and the plate allows a vertical adjustment adapted to the size of the animal.

The mask outlet can be connected with the gas evacuation apparatus to evacuate the waste gas and protect the safety of the operators and the environment.

3 different versions exist:

- a mask for mice adjustable in height and rotatable by 70 ° vertically and horizontally
- a mask for rats adjustable only in height
- a mask for rats adjustable in height and rotatable by 70 ° vertically and horizontally

Product features:**Available models**

Animal: Mouse (68057), Rat - up and down adjustment (68078), Rat - adjustment at different angles (68053)

