

## Endotracheal tube kit



for use with the R415 ventilator - up to 1kg small animal

**Product codes:**

Reference: 415K

### Product description:

This kit for use with the R415 ventilator includes below items:

- Endotracheal Tube for Rat, 14G, O.D.: 1.85mm
- Endotracheal Tube for Mice, 20G, O.D.: 0.9mm
- Luer Connector - Male, 3,2 mm
- Y-three way Valve - PP- 3,2 mm
- Silicone Tubing - ID: 3.0, OD: 6.0mm
- Reduced Connector - 3,2 to 6,3mm
- Silicone Tubing - ID: 6.0, OD: 9.0mm

### Product features:

### Available models