

## Mobile anesthesia machine



Compact and mobile solution for your anesthesia needs

**Product codes:**

Reference: R530

**Product description:**

Complete system, mounted on stand with handle, composed of a vaporizer, a flowmeter, inlet/outlet end caps and tube accessory kit.

This anesthesia system is mounted on a trolley equipped with 4 shelves for additional accessories sold separately: induction box, air pump, aspiration or weighing system for anesthetic gas cartridges ...

You can use this versatile system on rats, mice, hamsters, guinea pigs, rabbits, cats and other animals within 7kg.

You can choose your air supply (oxygen or air) and your anaesthetic product (isoflurane or sevoflurane).

**Product features:**

Weight: Less than 7 kg

Dimensions: 1335 mm H - 611 mm D

Temperature range: 10-35°C

Animal: Rats, mice, hamsters, guinea pigs, rabbits, cats and other animals within 7kg

Flow rate: 0.1-4L/min (oxygen flowmeter)

### **Available models**

Anaesthetic: Isoflurane (R530IE), Sevoflurane (R530SE)

