

## Accessories for RapID Tags

Applicator or Removal Tool for tag Rapid Lab



### Product codes:

Reference: BIO-RAPID

### Product description:

Applicator or Removal Tool for tag Rapid Lab

Applicator : Stainless Steel and Acetal , all components are autoclave safe at 121 degrees Celsius

Remove Tool : is not intended for autoclave but may be disinfected by alcohol

### Product features:

### Available models

RapId: removal tool (BIO-RAPIDREMOV), applicator (BIO-RAPIDAPPLI)