

# Peristaltic pump



Versatile and economical system

**Product codes:**

Reference: BS-900

## Product description:

This pump can be used as a **stand-alone pump** or as a **part of a network**. It is fully programmable and can be operated manually or with a PC through an RS-232 communication port.

2 versions are possible with different heads:

- **BS-9000G** has a 4 roller pump head that is designed for more accurate dispensing and higher pressure. It reaches a maximum of **75.19 ml/min**
- **BS-9004C** has a 2-roller pump head and a different shape from the 9000G pump head. It provides much higher flow rates up to **1038 ml/min**

## Product features:

Dimensions: 20 x 13.5 x 15 cm

Dimensions - Inner: Tubing of 1.6mm (1/16") for BS-9000G - 6.35 mm (1/4") for BS-9004C

Flow rate: 0.004 to 75.19 ml/min for BS-9000G - 0.047 to 1038 ml/min for BS-9004C

### Available models

Maximal flow rate: 75,19 ml/min (BS-900-0G), 1038 ml/min (BS-900-4C), Head kit : 2 Rollers, 3/16"i.d. x 5/16" (BS-900 -YB2), Head kit : 4 Rollers 1/16"i.d. x 3/16"o.d (BS-900 -YG4)

