

# Concentric tubing masks for mice and rat

Concentric Tubing Masks for mice and rats  
FOR



**Product codes:**

Reference: R510-M

## Product description:

Using Concentric Tubing design, sturdy and durable, compatibility between rat and mice achieved by replacing the inner core. The entire mask (including tubing) is made of non-metallic material suitable for MRI environments.

R510-M use together with R546Pro for positive evacuation

Included pipe

The mask can be fixed by fixed bracket (68625) or platform (R510-PGL).

## Product features:

## Available models