

# Ear Puncher

Ear punch with slanted edge, for identification and/or sample collection

**Product codes:**

Reference: EP-A7020

**Product description:**

This ear punch is used for identification and/or sample collection.

The angled edge allows cutting rather than puncturing, to minimize animal discomfort.

Requires very gentle pressure to cut the hole.

Length : 11cm

Biopsy sample 2mm

**Autoclavable** but must be lubricated after use for long and trouble-free use.

Made of stainless steel

**Product features:**

Diameter: Hole : 2mm

### Available models

