

# Red induction chamber for mice and rats

Red chamber, compatible with our anesthesia stations



## Product codes:

Reference: V105

## Product description:

The induction chamber is the best choice for rapid induction of anesthesia, by inhalation of anesthetic gases, before animal surgery. It takes only 2 to 5 minutes. It is also very useful when you can't use the masks with awake animals.

The anesthetic gas is inducted into the chamber, and exhausted to the connected charcoal canister (gas filter canister not included)

Dim : 31cm\*15cm\*19cm

## Product features:

Material: Plexiglass

Diameter: 6 mm input ID - 22 mm output ID

Thickness: 1.2 cm (upper cover) - 0.5 cm (other)

## Available models

