

Compact anesthesia machine



Compact solution for your anesthesia needs

Product codes:

Reference: TAIJI

Product description:

Complete system, mounted on stand with handle, composed of a vaporizer, a flowmeter, inlet/outlet end caps and tube accessory kit.

You can use this versatile system on rats, mice, hamsters, guinea pigs, rabbits, cats and other animals within 7kg.

You can choose your air supply (oxygen concentrator or air pump) and your anaesthetic product (isoflurane or sevoflurane).

You can have one animal in the chamber and another in the mask.

Product features:

Weight: Less than 7.5kg

Dimensions: 326mm (L) × 180mm (W) × 284mm (H)

Temperature range: 10-35°C

Animal: Rats, mice, hamsters, guinea pigs, rabbits, cats and other animals within 7kg

Flow rate: 0.1-4L/min (oxygen flowmeter)

Available models

Anaesthetic: Isoflurane (TAIJI-IE)

