

Thermo Controlled Surgery Platform



Thermo Controlled Surgery Platform with non-skid feet

Product codes:

Reference: ASS7T

Product description:

Thermo Controlled Surgery Platform -

Easily maintain animal body temperature during surgery!

Two flexible LED lights provide direct illumination with no shadows and no heat.

Thermo reflector and turbo heat sink to maximize heat transfer. Flex arm, non-obstructive rimless magnifying glass, 2.3x main lens, 4.0x inset lens.

Allows continuous eye contact with object when moving into magnification area.

Includes Limb restraints, 4 Space Gels for continuous heat. 2 Space Blankets to retain heat and non-skid feet for added stability.

Universal input power supply.

Size: 13.5" L x 9" W

Product features:

Available models

